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
WWW.COMPRABASCULA.COM

BALANZAS MULTIFUNCIONALES
BALANCES MULTIFONCTIONNELS
MULTI-FUNCTIONAL BALANCES

Manual de usuario | Manuel d'usager | User's guide

ANF-AND-ANG

El fabricante se reserva el derecho de modificar sin previo aviso las características de sus productos para introducir mejoras técnicas o cumplir con nuevas regulaciones oficiales. / Le constructeur se réserve le droit de modifier les caractéristiques de ses produits en vue d'y apporter des améliorations techniques ou de respecter de nouvelles réglementations. / The manufacturer reserves the right to modify the specifications of its products in order to make technical improvements or comply with new regulations.

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2. INTENDED USE

The balance you have acquired serves to determine the weighing value of the material to be weighed. It is intended to be used as a non-automatic balance, i.e. the material to be weighed is manually and carefully placed in the centre of the weighing plate. The weighing value can be read off after a stable weighing value has been obtained.

2.1. INAPPROPRIATE USE

Do not use the balance for dynamic weighing. In the event that small quantities are removed or added to the material to be weighed, incorrect weighing results can be displayed due to the stability compensation in the balance. (Example: slowly draining fluids from a container on the balance)

Do not leave a permanent load on the weighing plate. This can damage the measuring equipment.

Be sure to avoid impact shock and overloading the balance in excess of the prescribed maximum load rating (max.), minus any possible tare weight that is already present. This could cause damage to the balance.

Never operate the balance in hazardous locations. The series design is not explosion-proof.

Structural alterations may not be made to the balance. This can lead to incorrect weighing results, faults concerning safety regulations as well as to destruction of the balance.

The balance may only be used in compliance with the described guidelines. Varying areas of application/ planned use must be approved in writing.

2.2. BEFORE USE

Please read all instructions to familiarize yourself with the specification and operations before using and programming the scales.

2.3. ACCESSORIES

Adapter, Operating Manual.

3. INSTALLATION

The balances are designed in a way that reliable weighing results are achieved in common conditions of use.

You will work accurately and fast, if you select the right location for your balance.

3.1. INSTALLATION TIPS

- > Place the balance on a firm, level surface.
- > Avoid extreme heat and temperature fluctuation such as installation in the direct sunlight.
- > Protect the balance against direct draughts due to open windows and doors.
- > Avoid vibration during weighing.
- > Protect the balance against high humidity, moisture and dust.
- > Do not expose the device to extreme humidity for longer periods of time. Nonpermitted condensation (condensation of air humidity on the appliance) may occur if a cold appliance is taken to a considerably warmer environment. In this case, acclimatize the disconnected appliance for approx. 2 hours at room temperature after it has been disconnected from the mains.
- > Avoid static charging of the material to be weighed.
- > Major display deviations (incorrect weighing results) are possible if electromagnetic fields occur as well as due to static charging and instable power supply. It is then necessary to change the location.
- > In order to obtain exact results with the electronic balances, your balance must have reached the operating temperature. During this warming up time the balance must be connected to the power supply (mains, battery, rechargeable battery).

3.2. CALIBRATION

As the gravity value is not the same at every location on earth. To achieve best accuracy, the balance is highly recommended to be calibrated in compliance with the gravity at the place of use. This calibration process must be carried out for the first time of use, after each change of location as well as in case of environment temperature change.

To receive accurate measuring values it is also recommended to calibrate the balance periodically in weighing operation.

The calibration should be made with the recommended calibration weight (see specification).

3.3. PROCEDURE FOR CALIBRATION

Observe stable environmental conditions. A warming up time (see specification) is required for stabilization.

- > Turn on the balance by pressing the **ON/OFF** key.
- > Press **SET/TARE** key and keep it pressed, until [$\mathcal{L}^{\mathcal{R}}$] appears in the display for few seconds. After that the adjustment weight (see specification) appears and flashes in the display.
- > Place the adjusting weight in the centre of the weighing plate.
- > Now press the **SET/TARE** key. Short time later [$\mathcal{L}^{\mathcal{R}} F$] appears, then returns automatically to the normal weighing mode. Readings of the adjustment weight appear in the display.

An error during adjustment or the use of an incorrect adjusting weight will result in an error message [$\mathcal{L}^{\mathcal{R}} \mathcal{E}$]. Repeat adjustment.

Keep the adjustment close to the balance. Daily control of the weighing exactness is recommended for quality-relevant applications.

3.4. OPERATION USING A RECHARGEABLE BATTERY (OPTIONAL)

Lift-off the battery cover on the lower side of the balance. Connect the battery and then cover the battery door.

For battery operation the balance has an automatic switch-off function which can be activated or deactivated in the menu.

Procedures as follows:

- > Turn on balance using the **ON/OFF** key and wait for the zero display.
- > Press and hold the **FUNCTION** key until [$\overline{0000}$] appears on the display.
- > Press the **UNIT/COUNT** key or **HOLD/PRINT** key repeatedly until [\overline{OFF}] appears.
- > Confirm by pressing the **SET/TARE** key.

Use the **UNIT/COUNT** key and **HOLD/PRINT** key to choose between the following settings:

1. [$\overline{OFF ON}$] To save battery power, the balance switches off automatically if 3 minutes without weighing.
2. [$\overline{OFF OFF}$] Switch-off function deactivated.

Use the **SET/TARE** key to confirm selected setting. The balance returns to weighing mode.

If the batteries are exhausted, \overline{LO} appears in the display. Press **ON/OFF** key and replace the batteries immediately.

If the balance is not used for a longer time, take out the batteries and store them separately. Leaking battery liquid could damage the balance.

If there exists an optional rechargeable battery, it has to be connected in the battery compartment via a separate plug-in socket. The mains adaptor delivered with the rechargeable battery must be used.

3.5. USING ADAPTOR

- > Before using the external mains adaptor, make sure that the printed voltage must comply with the local voltage.
- > Only use original mains adaptor. The use of other adaptor is subject to approval by Giropès.

3.6. CONNECTION OF PERIPHERAL DEVICES

Before connecting or disconnecting of additional devices (printer, PC) to the data interface, always disconnect the balance from the power supply.

With your balance, only use accessories and peripheral devices by Giropès, as they are ideally tuned to your balance.

4. OPERATIONS

KEY	FUNCTION
ON/OFF	Turn on/off power of the balance
FUNCTION	Call up menu (keep key pressed until UNIT appears)
SET/TARE	- Confirm settings in the menu. - Save and exit menu. - Tare weight.
UNIT/COUNT	- Select weighing unit. - Change settings in the menu. - Parts counting.
HOLD/PRINT	- Display hold. - Calculate weighing data via interface. - Change settings in the menu.
ZERO	Zero balance.

4.1. WEIGHING

- > Turn on the balance using the **ON/OFF** key.
- > The balance will perform self-test for few seconds. When the zero display appears, the balance is ready for operations.

Important: If the display not showing zero, press the **ZERO** key.

- > Place objects onto the center of pan. Take care that the weighed material should not touch the balance housing.
- > Now the weight is displayed. Once the reading is stable, the weighing unit (e.g. *g* or *kg*) will appear on the right-hand side of the display.
- > If the objects are heavier than the weighing capacity, the display will show *Error* (=overload).

4.2. TARING

- > Turn on the balance using the **ON/OFF** key and wait for the *ZERO* display.
- > Put the tare vessel on the weighing plate and press the **SET/TARE** key. The balance display goes to zero. The weight of the vessel is now saved.
- > The taring process can be repeated any number of times, e.g. when adding several components for a mixture (adding).
- > The limit is reached when the total weight reaching full capacity of the balance. After removing the taring vessel, the total weight is displayed as negative reading.
- > After finishing the weighing process, press the **SET/TARE** key again, the balance display goes back to *ZERO*.

4.3. PRE-TARE FUNCTION

Using this function the weight of a tare vessel can be stored.

This value also remains saved if the balance meanwhile has been switched off and switched on again.

- > Turn on balance using the **ON/OFF** key and wait for the zero display.
- > Put tare vessel on the weighing plate. Press the **UNIT/COUNT** key until balance displays *5 PCS*.
- > Press the **UNIT/COUNT** key 6 times until *PRE* flashes on the display. By actuating the **SET/TARE** key, the current weight on the weighing plate is saved as PRE-Tara weight.
- > To switch off this function, unload the weighing plate. Press **ZERO** key if display is not showing zero.
- > Press the **UNIT/COUNT** key until balance displays *5 PCS*.
- > Press the **UNIT/COUNT** key 6 times until *PRE* flashes on the display. By actuating the **SET/TARE** key, the stored PRE-Tare weight is deleted.

4.4. PLUS/MINUS WEIGHING

This function can be used for unit weight control, fabrication control etc.

- > Turn on balance using the **ON/OFF** key and wait for the *ZERO* display.
- > Put the nominal weight on the weighing plate and tare to zero by pressing the **ZERO** key. Remove the nominal weight.
- > Put the test objects subsequently on the weighing plate, the respective deviation from the nominal weight is displayed with the sign “+” o “-”.
- > According to the same procedure also packages with the same weight can be produced, referring to a nominal weight.
- > Back to weighing mode by pressing the **ZERO** key.

4.5. PARTS COUNTING

- > Turn on balance using the **ON/OFF** key and wait for the zero display.
- > Press the **UNIT/COUNT** key until balance displays *5 PLS*. By pressing the **UNIT/COUNT** key several times, more reference quantities 10, 20, 25 and 50 can be called up.
- > Place the same number of pieces to count on the weighing plate as the set reference quantity requires.
- > Confirm with **SET/TARE** key.
- > The balance is now in parts counting mode counting all units on the weighing plate.
- > By pressing the **UNIT/COUNT** key the balance returns to the weighing mode and displays the weight of the counted units.

Important: The larger the reference quantity, the more accurate the parts counting.

*If the weight is less than the minimum unit weight, the display appears *Er 1*. Use **UNIT/COUNT** key to return to weighing mode. For minimum unit weight, see table “Specifications”.*

Tare vessels can also be used for piece counting. Pre-tare the weight of the tare vessel prior to piece counting.

4.6. NET-TOTAL WEIGHING

It is useful if a mixture of several components is weighed into a tare vessel and finally the sum weight of all weighed components is necessary for control purposes (net-total, i.e. the weight of the tare vessel).

EXAMPLE

- > Turn on balance using the **ON/OFF** key and wait for the *ZERO* display.
- > Put tare vessel onto weighing plate, tare with **ZERO** key to zero.
- > Weigh component. Press the **SET/TARE** key. Display shows zero with the memory icon, triangle, on the left side of the display.
- > Weigh component. Press the **SET/TARE** key. Net-total, sum weight of the components and will be displayed.
- > Tare to zero using the **SET/TARE** key. Weigh component. Press the **SET/TARE** key. Net-total, sum weight of components and will be displayed.
- > If necessary, repeat above procedures to complete the formula up to the desired final value. Repeat.
- > Back to weighing mode by pressing the **SET/TARE** key.

4.7. PERCENTAGE WEIGHING

Display symbol: %

Percent weighing allows to display weight in percent, in relation to a reference weight.

- > Turn on balance using the **ON/OFF** key and wait for the zero display.
- > Press the **UNIT/COUNT** key until balance displays *5 PLS*. Pressing the **UNIT/COUNT** key several times, until *100%* is displayed.
- > Place the reference item on the weighing pan.
- > Press **SET/TARE** key, the weight of the item is taken over as reference (100%).
- > Now you can place the test objects onto the weighing plate; the percentage is displayed.
- > Back to weighing mode by pressing the **UNIT/COUNT** key.

4.8. WEIGHING UNITS

Turn on balance using the **ON/OFF** key and wait for the zero display.

Press and hold the **FUNCTION** key until *UNIT* appears on the display.

Press **SET/TARE** key, the selected unit is displayed.

Use the **UNIT/COUNT** key or **HOLD/PRINT** key to select between different units (see following table).

Press the **PROGRAMAR** key to select the default weighing unit. It will be displayed whenever the balance is turned on.

UNIT	DISPLAY	CONVERSION FACTOR 1g=
Gram	g	1
Ounce	oz	0.035273962
Troy ounce	ozt	0.032150747
Pound	lb	0.002204623
Tael Hongkong	tlh	0.02671725
Tael Taiwan	tlt	0.0266666
Grain (depends on model)	gn	15.43235835
Pennyweight (depends on model)	dwt	0.643014931
Momme	mom	0.2667
Tola	tol	0.085733338
Carat (depends on model)	ct	5
Free Factor *	FFA	xx.xx

* In order to enter an own conversion factor:

- > Press the **UNIT/COUNT** key or **HOLD/PRINT** key as explained above until *FFA* is displayed.
- > Press the **SET/TARE** key to confirm the selection. The last digit begins to flash.
- > Press the **UNIT/COUNT** key to increase the value; press the **HOLD/PRINT** key to decrease the value.
- > Press the **ZERO** key to set the next digit.
- > When all digits are set, press the **SET/TARE** key to save this value.
- > The “Free Factor” is now the current weighing unit.

4.9. DISPLAY BACK LIGHT

The display back light can be set according to the applications:

- > Turn on balance using the **ON/OFF** key and wait for the zero display.
- > Press and hold the **FUNCTION** key until *UNIT* appears on the display.
- > Press the **UNIT/COUNT** or **PRESIONAR/IMPRIMIR** key several times until *bl* is displayed. Press **SET/TARE** to enter back light setting.
- > Press the **UNIT/COUNT** or **PRESIONAR/IMPRIMIR** key to select setting based on below table.

DISPLAY	ADJUSTMENT	FUNCTION
<i>bl on</i>	Back Light ON	Good contrast display which can also be read in darkness
<i>bl off</i>	Back Light OFF	Battery saving
<i>bl ch</i>		Back light will be switched off automatically if no weight change for 10 seconds. Battery saving.

Press the **SET/TARE** key to confirm setting.

4.10. ANIMAL WEIGHING FUNCTION

The balance has an integrated animal weighing function (mean value calculation).

With this function it is possible to weigh domestic or small animals, even they do not stand still on the weighing plate.

Note: If the animal moves too much, an exact weighing will not be possible.

In the menu, the animal weighing function can be switched on or off. To achieve this, follow the sequence of operations below:

- > Turn on balance using the **ON/OFF** key and wait for the *ZERO* on display.
- > Press and hold the **FUNCTION** key until *UNIT* appears on the display.
- > Press the **UNIT/COUNT** or **PRESIONAR/IMPRIMIR** key several times until *ANL* is displayed.
- > Press **SET/TARE** key to display current setting.
- > Press the **UNIT/COUNT** or **PRESIONAR/IMPRIMIR** key to select setting based on below table.

DISPLAY	FUNCTION
<i>ANL OFF</i>	Animal weighing function is switched off.
<i>ANL 3</i>	Weighing value calculation for 3 seconds before displays weight.
<i>ANL 5</i>	Weighing value calculation for 5 seconds before displays weight.
<i>ANL 10</i>	Weighing value calculation for 10 seconds before displays weight.
<i>ANL 15</i>	Weighing value calculation for 15 seconds before displays weight.

Press the **SET/TARE** key to confirm setting.

OPERATION

Turn on balance using the **ON/OFF** key and wait for the zero display.

Put the weighing object (animal) on the weighing plate and press the **SET/TARE** key. In the display the preselected time is displayed in seconds and counting down. During this period, the balance is calculating the mean value. The calculated weight is then displayed and hold on the screen.

By pressing the **SET/TARE** key, the balance switches between normal weighing mode and animal weighing mode.

4.11. DISPLAY READING HOLD FUNCTION

The balance has the display reading hold function so that important reading can be captured for a period. Press the **HOLD/PRINT** key once during weighing, display reading can be frozen. Press the **HOLD/PRINT** key again to change back to normal weighing mode.

4.12. AUTO ZERO TRACKING FUNCTION

The auto zero tracking function is used to tare out small variations in weight automatically. In some applications, small quantities are removed or added to the material to be weighed, this feature can affect the weighing results due to the "stability compensation" function in the balance. (Example: Slowly draining fluids from a container on the balance).

When application involves small variations of weight, it is advisable to switch off this function. However, the weight reading will look less stable.

In weighing mode, press and hold the **FUNCTION** key until *UNIT* appears on the display.

Press the **UNIT/COUNT** or **HOLD/PRINT** key several times until *LR* is displayed.

Press **SET/TARE** key to display current setting.

Press the **UNIT/COUNT** or **HOLD/PRINT** key to select setting based on below table.

- > *LR* **on** Function activated
- > *LR* **off** Function deactivated

Press the **SET/TARE** key to confirm selection. The balance returns to weighing mode.

4.13. RESET TO FACTORY SETTING

This function resets all balance settings to factory setting.

In weighing mode, press and hold the **FUNCTION** key until *UNIT* appears.

Press the **UNIT/COUNT** key or **HOLD/PRINT** key several times until *rSL* is displayed.

Press the **SET/TARE** key to display current setting.

Select the desired settings by pressing the **UNIT/COUNT** key or **HOLD/PRINT**.

- > *rSL* **SI** Balance will be reset to factory setting
- > *rSL* **no** Balance keeps its individual setting

Press the **SET/TARE** key to confirm selection. The balance returns to weighing mode.

5. CAUTIONS

- > The balance is not a water proof unit. Splash of water into the scale may cause permanent damage to balance.
- > Do not over load the balance which may cause permanent damage of the sensor
- > Always use the balance with care since this is a very precision device.
- > Do not use balance for dynamic weighing. In the event that small quantities are removed or added to the material to be weighted, incorrect weighting results can be displayed due to the “stability compensation” in the balance.
- > Do not leave permanent load on the weighing plate. This may damage the measuring system.
- > Never operate balance in explosive environment. The version is not explosion protected.
- > The structure of the balance may not be modified. This may lead to incorrect weighing results, safety-related faults and destruction of the balance.

6. WARRANTY

This scale is warranted against defects of manufacturing and materials for a period of 1 year, from the delivery date.

During this period, GIROPÈS SL will take charge of repairing the scale.

This warranty does not cover defects or damaged caused by misuse, overloading or improper installation contrary to the recommendations described in this manual.

This warranty does not cover shipping costs for the reparation of the balance.

NOTAS - NOTES

NOTAS - NOTES

Baxtran

marca propiedad de /
marque propriété de /
trade mark propriety of:

Pol. Empordà Internacional

C/ Molló, 15-16
17469 VILAMALLA
(Girona) SPAIN

T. (34) 972 527 212
F. (34) 972 527 211